

WelchAllyn®



Welch Allyn
Connex® Spot Monitor

Welch Allyn Connex Spot Monitor

Accurate, automatic vitals sent to your EMR from the point of care to help enhance staff productivity, patient outcomes and risk reduction

The only monitor to turn to for truly wireless data transfer, accurate vital signs measurement across your patient population and compatibility with your EMR now and in the future.

Accurate readings across your patient population

- > Help ensure accurate hypertension diagnosis without slowing down your workflow with unique SureBP® technology and blood pressure averaging.
 - Research shows that more than half of patients could be misdiagnosed in a hypertensive category using the single-reading method.¹
- > Provide SpO₂ accuracy in challenging conditions including motion, thick fingers, dark skin tone and low perfusion with Masimo® SET® and Nonin® PureSAT® pulse oximetry technology.²

Improve data accuracy, efficiency and mobility

- > The Connex Spot Monitor is designed to help eliminate data entry errors in the patient chart with wireless transfer of vital signs.
- > Increase staff efficiency by eliminating typing into patient records.
- > Provide freedom from wires and cables, giving you the flexibility to easily work in limited space.



Configurable design

allows you to customize the monitor to the needs of your practice and your patients

The easy-to-use, vivid touchscreen display is perfect for the physician office environment. Select from these options to get exactly the monitor you need:



Thermometry

- SureTemp® Plus Thermometry
- Braun ThermoScan® PRO 6000 Ear Thermometry

Pulse Oximetry

- Masimo SET SpO₂
- Nonin PureSAT SpO₂
- Nellcor® OxiMax® SpO₂

Mounting and Mobility Options

- **Integrated Wall System** combines physical assessment and digital vital signs for a complete patient exam workflow
- **Accessory Power Management Stand** features an integrated battery to ensure the device gets you through a busy day
- **Standard Basket Stand** provides the simple portability you need with a minimal footprint



Connectivity for Ambulatory Care

- > Welch Allyn closely partners with more than 90 leading EMR providers—including Epic®, athenahealth® and NextGen®—to ensure reliable device connectivity and to help you leverage your EMR investment.
- > Connect your Connex Spot Monitor to scales from Detecto®, Health o meter® and SECA to accurately document height, weight and BMI.

Protect your investment with service programs

- > **Partners in Care™ Comprehensive Partnership Programs** help minimize the cost and complexity of servicing and maintaining your monitors.
- > **Enterprise Deployment Option** enables your monitors to arrive configured to your specification right out of the box.

¹Robert Smith, MD. [Blood pressure averaging methodology: Decreasing the rate of misdiagnosing hypertension](#). 2014.

²Nonin Medical, [Patient conditions can be unpredictable. Pulse oximeter accuracy shouldn't be](#). 2014.

Configurations Available

CONNEX SPOT MONITOR	7100	7300
Parameters		
SureBP Blood Pressure	●	●
Blood Pressure Averaging	●	●
SureTemp Plus Thermometry	○	○
Braun ThermoScan PRO 6000 Ear Thermometry	○	○
Masimo or Nellcor OxiMax SpO ₂		○
Nonin SpO ₂	○	○
Communications		
Wireless (Bluetooth®)		●
USB, Ethernet	●	●
Accessories		
Standard Basket Stand	○	○
Accessory Power Management Stand	○	○
Green Series™ 777 Integrated Wall System	○	○

● Included ○ Optional

ACCESSORIES

7000-APM	Connex Spot Accessory Power Management Stand
7000-MS3	Connex Spot Classic Mobile Stand
7000-GCX	Connex Spot Wall Channel
77794-MCSM	Green Series 777 Integrated Wall System with Coaxial Ophthalmoscope
77794-2MCSM	Green Series 777 Integrated Wall System with PanOptic Ophthalmoscope
7000-PS	Power Supply (required if not ordering a stand or wall mount)
BATT22	Connex Spot Monitor Lithium-Ion Battery
BATT99	Accessory Power Management Stand Lithium-Ion Battery
106275	Connectivity Accessory Kit USB Cable
4500-35	Blood Pressure Hose with FlexiPort® (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
02895-100	SureTemp Plus Thermometer Rectal Probe and Well
718584	Braun PRO 6000 Ear Thermometer Tether
3278-010	Nonin Reusable Adult Sensor (2 meters)
2360-010	Nonin Reusable Pediatric Sensor (2 meters)
S1-CSM-5	Connex Spot Monitor Comprehensive Partnership Program, 5 years*
S2-CSM-5	Connex Spot Monitor Biomed Partnership Program, 5 years*

*One- and two-year service plans also available.



Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM



Welch Allyn®