

# Video Laryngoscope

PRODUCT CODE: UB-VL3R



TEL: 0 100 350 238

FAX : 086 477 9020

E-MAIL: [info@ububelemedical.co.za](mailto:info@ububelemedical.co.za)

## FEATURES

- **Display size:** 3.5"LCD Full View
- **Resolution:** 640\*480 (RGB, 2M pixel high resolution)
- **Camera Function:** Get serial camera shooting by pressing the button on the handle for training and presenting
- **Uploading:** Upload picture and video by USB port
- **Illumination:** LED
- **Image Output:** Connect to external monitor and convenient for training and presentation.
- **Temperature(Working condition):** 5°C~ + 40°C
- **AP(Working condition):** 700hPa to 1060hpa
- **Unique Anti-fog Function:** Get anti-fog function upon powering without preheating
- **Handle:** Comfortable ergonomic design ; anti-microbial
- **Blade:** Reusable. blade can be sterilized by cleaner and reused for one thousand times; Disposable one is only for one time and then thrown away





# Video Laryngoscope

## SPECIFICATIONS

- **N.W(g/pc):** 350g
- **G.W(g/pc):** 1.5KG
- **Package Dimensions:** 27\*27.5\*12CM
- **Battery Type:** Rechargeable Lithium Battery , :3 years+
- **Long working time:** 4 hrs continuous working time
- **Voltage:** 3.7V
- **Capacity:** 3200mAh
- **Charging times &time:** ≥300 times ;<8H
- **Charging port(Power Adapter):** Micro USB
- **Input::** 100-250V,50Hz
- **Output:** 5V,1000mA
- **Temperature:** "-10°C~ + 45°C"
- **Humidity:** ≤93%
- **AP:** 500hPa to 1060hPa



TEL: 0 100 350 238

FAX : 086 477 9020

E-MAIL: [info@ububelemedical.co.za](mailto:info@ububelemedical.co.za)