

EARIGATOR™ EAR WASH SYSTEM

earigator™



Designed by an Otologist, the Earigator™ combines the functions of an otoscope and irrigation into one, providing the most advanced means of cerumen control.

Fully automated temperature and pressure controls make the Earigator™ a safe and highly effective ear wax removal system.

Features

- Fully automatic temperature and pressure controls
- Non-physician, one person operation
- 3-5 minute procedure with maximum safety
- LED illumination with magnification
- System always on standby
- Highly effective and pain free



TEL: 0 100 350 238

FAX : 086 477 9020

E-MAIL: info@ububelemedical.co.za

EARIGATOR™ EAR WASH SYSTEM

Specifications

Operating Temperature	37 ± 2°C (98 ± 3°F)
Operating Pressure	0.8 ± 0.1 bar (12 ± 2 psi)
Water Flow Rate	1.25 L/min (0.33 gal/min)
Warm-Up Time	5-10 minutes
Water Capacity	4 liters (1 gallon)
Electrical	100-120 VAC, 50/60 Hz, 10 A (model NT0001-1) 220-240 VAC, 50/60 Hz, 5 A (model NT0001-2)
Dimensions	32L x 27W x 40H cm (12L x 10W x 16H in)

Base Unit - pressure, temperature & flow controls, water reservoir, heater, pump, and digital display

Hand Piece - flow trigger, water nozzle, magnifying glass, and LED lights

Accessories - catch basin, disposable nozzles, and stainless steel rolling cart



Ordering Information

- **UB-NT0001-2** Earigator™ 220V Unit
- **UB-MI0001** Optional Cart
- **UB-CA0005** Catch Basin
- **UB-FO0002-2** Disposable Nozzles (100 PC)



TEL: 0 100 350 238

FAX : 086 477 9020

E-MAIL: info@ububelemedical.co.za