

MCW Health and Fitness Scale



MCW 300L



Technical Specifications

Model	MCW 300L (Digital)
Maximum capacity	300kg
Readability	0.05kg
Repeatability (S.D.)	0.1kg
Linearity (±)	0.1kg

Applications

- Weighing


Features MCW 300L (Digital)

- Large, clear LED display
- Colour coded, sealed keypad
- Hold function
- Swivel foot and armrests
- Swivel wheels for easy movement and foot brake
- Full range tare
- Weighing in kg and lb
- RS-232 Interface
- Adapter / Rechargeable Battery Power
- BMI - Body Mass Index

General Specifications

Model	MCW 300L (Digital)
Units of measurement	Kilogram (kg), Pound (lb)
Calibration	External
Display	LED display
Interface	RS-232
Body mass index	Standard
Power supply	12vDC 500mA adapter / Internal Rechargeable Battery
Operating Temperature	5° to 40°C
Housing	Base: mild steel Indicator: ABS plastic Chair: plastic
Overall dimensions (wxdxh)	675mm x 875mm x 950mm
Net weight	25kg



 010 035 0238
 086 477 9020
 info@ububelemedical.co.za
 www.ububelemedical.co.za