



FEATURES

- BMI (body mass index) function available on the MDW 300L
- The MDW 300L offers RS-232 interface for fast communication with printers and computers
- MDW 250L and 300L feature a hold function to retain the weight on the display
- Use the MDW 250L and 300L with the included AC adapter or rechargeable battery
- Height measurements from 60cm to 212cm
- Anti-slip platform ensures safety stepping on and off the scale
- Integrated wheels permit easy movement (MDW300L only)
- Height rod facilitates measuring the height of children and adults while weighing
- Tare bar included with MDW200M only

Affordable, easy-to-operate scale measures height and weight

GENERAL SPECIFICATIONS

- Weighing units: kg, lb (L models only)
- Display: Graduated beam (M models only)
- Calibration: External calibration
- Operating Temperature: 5° to 40°C
- Power Supply: 15VDC 50/60Hz 500mA adapter / Internal rechargeable battery (L models only)
- Display: LED with 22mm-high digits (L models only)
- Interface: RS-232 (MDW300L only)

APPLICATIONS

- Weighing

WEIGHT AND DIMENSIONS

- **MDW 160M:** Net Weight - 15kg, Overall Dim. - 530x275x1485mm (wxdxh)
- **MDW 200M:** Net Weight - 15kg, Overall Dim. - 530x275x1485mm (wxdxh)
- **MDW 250L:** Net Weight - 15kg, Overall Dim. - 530x275x1485mm (wxdxh)
- **MDW 300L:** Net Weight - 13kg, Overall Dim. - 570x275x1310mm (wxdxh)

| Model | Capacity | Readability | Pan Size |
|-----------------|----------|-------------|-----------|
| MDW 160M | 160kg | 0.1kg | 375x275mm |
| MDW 200M | 200kg | 0.1kg | 375x275mm |
| MDW 250L | 250kg | 0.1kg | 375x275mm |
| MDW 300L | 300kg | 0.05kg | 375x275mm |



ACCESSORIES for MDW300L only

| Item number | Description |
|-------------|---|
| 1120011156 | ATP thermal printer |
| 3126011263 | ATP thermal printer paper |
| 3126011281 | ATP thermal printer paper (pack of 10) |
| 3014011014 | RS-232 cable |
| 3074010507 | RS-232 to USB adapter (must order corresponding RS-232 cable) |
| 600002028 | Adam DU data collection program |
| 700660290 | Calibration certificate |



010 035 0238
 086 477 9020
 info@ububelemedical.co.za
 www.ububelemedical.co.za