

MDW Health and Fitness Scales



Technical Specifications

Model	MDW 250L (Digital)	MDW 300L (Digital)
Maximum capacity	250kg	300kg
Readability	0.1kg	0.05kg
Repeatability (S.D.)	0.2kg	0.1kg
Linearity (±)	0.2kg	0.1kg
Pan size (w x d)	275mm x 375mm	

MDW 250L

Applications

- Weighing
- Measuring height



MDW 300L

Features MDW 250L & 300L (Digital)

- Anti-slip platform
- Large, clear LED display
- Sealed Keypad
- Hold Function
- Adapter / Rechargeable Battery Power
- Robust painted mild steel base
- Weighing in kg and lb
- Adjustable height rod (measures up to 210cm)
- Enhanced height rod (MDW 300L model only)
- BMI-Body Mass index (MDW 300L model only)
- Wheels for easy moving (MDW 300L model only)
- RS-232 Interface (MDW 300L model only)

General Specifications

Model	MDW 250L (Digital)	MDW 300L (Digital)
Units of measurement	Kilogram (kg), Pound (lb)	Kilogram (kg), Pound (lb)
Calibration	External	
Display	LED display	
Interface	-	RS-232 MDW 300L only
Body mass index (BMI)	-	MDW 300L only
Power supply	12vDC 500mA adapter / Internal Rechargeable Battery	
Operating Temperature	5° to 40°C	
Housing	Base: mild steel with anti slip platform Indicator: ABS plastic	
Overall dimensions (wx dxh)	535mm x 275mm x 940mm	570mm x 275mm x 1310mm
Net weight	13kg	



 010 035 0238
 086 477 9020
 info@ububelemedical.co.za
 www.ububelemedical.co.za