



MTB Baby and Toddler Scale

Lightweight and portable, the MTB baby and toddler scale offers excellent value and versatility for pediatricians and other healthcare providers. With a capacity of 20kg, this scale boasts the flexibility to weigh newborn babies to toddlers on the same device. The MTB features quick start-up and easy operation, so minimal training is needed. The cleverly designed ergonomic cradle provides comfortable support for babies and locks securely in place. The cradle removes easily to reveal a non-slip rubber mat underneath, simplifying the task of weighing toddlers. A memory function retains the most recent weight taken, allowing the user to record results after safely removing the child from the scale.



MTB Baby and Toddler Scale

Features

- High-contrast display is easily visible from a distance
- Memory features can recall the last weight used allowing the user to focus on the baby or toddler
- Cradle can be removed for weighing toddlers on the non-slip rubber mat
- Ergonomic cradle provides comfortable support
- Lightweight design simplifies transport and portability
- Auto power-off to save energy
- External calibration allows for verification and adjustment with weights
- Simple operation allows quick startup and use with little or no training
- Operates with AC adapter (included) or batteries (not included)

Applications

- Weighing

Model	MTB20
Capacity	20kg
Readability	5g
Repeatability (S.D.)	5g
Linearity (+/-)	10g
Pan Size	285x560mm
Weighing Units	kg, lb, lb:oz
Stabilisation Time (sec)	2
Calibration	External calibration
Display	LCD with 29mm-high digits
Power Supply	6VDC 500mA Adapter
Operating Temperature	10° to 40°C
Housing	ABS
Overall Dim.	560x350x140mm (wxdxh)
Net Weight	3.2kg



010 035 0238

086 477 9020

info@ububelemedical.co.za

www.ububelemedical.co.za